



Essential Medical Products

LTC Bariatric & Ultra Bariatric Mattress

• **Bariatric mattress:** Dual-layer construction designed for heavier clients up to 500 lbs.

• **Ultra Bariatric mattress:** Tri-layer construction accommodates clients up to 700 lbs.



Patient Considerations for this Mattress:

- Pressure Distribution** › Good choice for individuals at risk for pressure sore development
- Firmness** › Recommended for individuals requiring a medium-firm mattress
- Temperature** › Functional choice for prevention of heat build-up and sweating
- Transfers** › No transfer borders - reduced stability for unassisted transfers
- Bed Mobility** › Good choice for ease of movement on bed surface
- Other** › Best choice for bariatric individuals

Specifications:

- Top Foam Layer** › 3" medium density foam comfort layer
- Core Foam** › Bariatric has a multi-zone firm core / Ultra Bariatric has a dual-layer firm core
- Transfer Borders** › No exterior transfer borders
- Cover** › Fluid-proof, non-allergenic, antibacterial, antimicrobial cover
Also available with organic cotton cover (not fluid-proof)
- Misc** › Designed for use on all hospital beds

Standard Sizes: (Custom sizes also available by quotation.)

- 80" x 35" x 7"
- 80" x 36" x 7"
- 80" x 39" x 7"
- 80" x 48" x 7"
- 80" x 54" x 7"
- 80" x 59" x 7"

Recommended Weight Capacity:

Bariatric Mattress - 500 lb / Ultra Bariatric Mattress - 700 lb

100% Customer Satisfaction Guarantee:

The warranty period for this product is 1 year.

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.

