


## LTC 3500 Active Care Mattress

Designed for clients with muscle weakness who find it difficult to move in bed. Combines a top contour cut comfort layer to reduce shear and multi-zoned core enhancing pressure distribution over vulnerable areas.



### Patient Considerations for this Mattress:

- |                              |   |
|------------------------------|---|
| <b>Pressure Distribution</b> | › Functional choice for individuals at risk for pressure sore development |
| <b>Firmness</b>              | › Recommended for individuals requiring a medium-firm mattress            |
| <b>Temperature</b>           | › Functional choice for prevention of heat build-up and sweating          |
| <b>Transfers</b>             | › Good choice for bedside stability, preventing entrapment                |
| <b>Bed Mobility</b>          | › Best choice for ease of movement on bed surface                         |

### Specifications:

- |                         |   |
|-------------------------|---|
| <b>Top Foam Layer</b>   | › 2" plush foam comfort layer   |
| <b>Core Foam</b>        | › Multi-zone medium-firm core   |
| <b>Transfer Borders</b> | › 4" high exterior foam adds firmness along edges   |
| <b>Cover</b>            | › Fluid-proof, non-allergenic, antibacterial, antimicrobial cover<br>Also available with organic cotton cover (not fluid-proof) |
| <b>Misc</b>             | › Designed for use on all hospital beds.  |

### Standard Sizes: (Custom sizes also available by quotation.)

- 80" x 35" x 6"
- 80" x 36" x 6"
- 80" x 39" x 6"
- 80" x 48" x 6"
- 80" x 54" x 6"
- 80" x 59" x 6"

### Recommended Weight Capacity:

350 lbs

### 100% Customer Satisfaction Guarantee:

The warranty period is five years; two years full plus three years pro-rated.

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.