



LTC 4000 Ultra 55 Mattress

Ultra Stable Firm Memory Foam

Firm viscoelastic foam plus a high resiliency transitional layer with full height transfer borders ensure stability during transfers and when performing ADL's. Helps prevent entrapment between mattress and bedside rails.



Patient Considerations for this Mattress:

- Pressure Distribution** › Best choice for individuals at risk for pressure sore development
- Firmness** › Recommended for individuals requiring a firm mattress
- Temperature** › Good choice for prevention of heat build-up and sweating
- Transfers** › Best choice for bedside stability, preventing entrapment
- Bed Mobility** › Viscoelastic top layer may impair bed mobility for some individuals
- Other** › Good choice for heavier individuals or individuals preferring a firm mattress

Specifications:

- Top Foam Layer** › 2" medium-plush viscoelastic foam comfort layer plus 2" high resiliency transitional layer
- Core Foam** › Firm core
- Transfer Borders** › Full height exterior foam adds firmness along edges
- Cover** › Fluid-proof, non-allergenic, antibacterial, antimicrobial cover
Also available with organic cotton cover (not fluid-proof)
- Misc** › Designed for use on all hospital beds.

Standard Sizes: (Custom sizes also available by quotation.)

- 80" x 35" x 6"
- 80" x 36" x 6"
- 80" x 39" x 6"
- 80" x 48" x 6"
- 80" x 54" x 6"
- 80" x 59" x 6"

Recommended Weight Capacity:

450 lbs

100% Customer Satisfaction Guarantee:

The warranty period is five years; two years full plus three years pro-rated.

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.