


## LTC 8000 Gel-In Mattress

For Heat Regulation and Energy Absorption

Advanced hybrid gel/foam technology with cooling gel infused into the viscoelastic foam top layer which results in a perfect balance of heat regulation and energy absorption. Maximum cooling ability for clients who sweat due to overheating.



### Patient Considerations for this Mattress:

- |                              |   |
|------------------------------|---|
| <b>Pressure Distribution</b> | › Best choice for individuals at risk for pressure sore development   |
| <b>Firmness</b>              | › Recommended for individuals requiring medium firmness               |
| <b>Temperature</b>           | › Best choice for prevention of heat build-up and sweating            |
| <b>Transfers</b>             | › Best choice for bedside stability and preventing entrapment         |
| <b>Bed Mobility</b>          | › Viscoelastic top layer may impair bed mobility for some individuals |

### Specifications:

- |                         |   |
|-------------------------|---|
| <b>Top Foam Layer</b>   | › 2" gel infused viscoelastic foam comfort layer  |
| <b>Core Foam</b>        | › Multi-zone medium-firm core   |
| <b>Transfer Borders</b> | › 6" high exterior foam adds firmness along edges   |
| <b>Cover</b>            | › Fluid-proof, non-allergenic, antibacterial, antimicrobial cover<br>Also available with organic cotton cover (not fluid-proof) |
| <b>Misc</b>             | › Designed for use on all hospital beds   |

### Standard Sizes: (Custom sizes also available by quotation.)

- |                  |                  |
|------------------|------------------|
| • 80" x 35" x 6" | • 80" x 48" x 6" |
| • 80" x 36" x 6" | • 80" x 54" x 6" |
| • 80" x 39" x 6" | • 80" x 59" x 6" |

### Recommended Weight Capacity:

350 lbs

### 100% Customer Satisfaction Guarantee:

The warranty period is five years; two years full plus three years pro-rated.

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.