




Essential Medical Products

LTC 9200 T-Style Configurable Mattress

Pressure Re-distribution Focused on Coccyx and Heels



Unique design configuration allows for optimal pressure re-distribution focusing on high at-risk areas while reducing the number of air/gel inserts needed.

This mattress when used in conjunction with either air or gel inserts creates an optimal therapeutic environment for skin healing and protection.

Patient Considerations for this Mattress:

- Pressure Distribution** † Best choice for individuals at risk for pressure sore development at the coccyx and heels
- Firmness** † Top comfort layer of medium firm viscoelastic foam in non high-risk areas
- Temperature** † Best choice for prevention of heat build-up when air or gel inserts used
- Transfers** † Best choice for bedside stability and prevention of entrapment with full height transfer borders
- Bed Mobility** † Viscoelastic foam top layer may impair bed mobility for some individuals

Specifications:

- Top Foam Layer** † Top surface comprised of medium-firm viscoelastic foam
- Core Foam** † Medium-firm polyurethane foam
- Transfer Borders** † Transfer borders add firmness along sides
- Cover** † Fluid –proof, non-allergenic, antibacterial, and anti-microbial (An additional all-way stretch Emersion “shower-cap” style cover is also available)
- Misc** † Designed for use on all hospital beds

Standard Sizes: (Custom sizes also available by quotation.)

- 80" x 35" x 6"
- 80" x 36" x 6"

Recommended Weight Capacity:

350 lbs

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.

100% Customer Satisfaction Guarantee:

The warranty period for this product is 5 years.

