

# LTC 9000 Progressive Care Mattress

## Adjustable Pressure Redistribution

Versatile design allows multiple configurations to meet specific patient needs. Use a combination of interchangeable foam and/or air inserts to optimize care and comfort. A removable 'single envelope' cover protects internal components, ensuring easy cleaning.

Side View



### Patient Considerations for this Mattress:

- Pressure Distribution** › Best choice for individuals at risk for pressure sore development when air-filled or gel inserts used in high risk areas that may change location over time
- Firmness** › Firmness can be adjusted to suit individual need
- Temperature** › Best choice for prevention of heat build-up and sweating when used with add gel or air-filled inserts
- Transfers** › Best choice for bedside stability and preventing entrapment when used with full height transfer borders
- Bed Mobility** › Viscoelastic top layer may impair bed mobility for some individuals

### Specifications:

- Top Layer Inserts** › Choose a combination of:
  - a) Gel-in foam - Gel infused viscoelastic foam provides superior pressure distribution plus the cooling affect of gel
  - b) Medium-firm viscoelastic foam - conforms to body contours, providing superior pressure redistribution
  - c) Medium-soft polyurethane foam - provides firm support
  - d) Both air-filled and gel inserts are available for your dealer  
Contact your local supplier or EMP for details
- Core Foam** › Medium-firm core
- Transfer Borders** › Choose 3" or 6" height. 3" transfer border allows top foam layer to extend over full mattress width. 6" transfer border adds firmness along sides
- Cover** › Fluid-proof, Non-allergenic, antibacterial, anti-microbial single envelope cover with three sided double slider zipper
- Misc** › Designed for use on all hospital beds

### Standard Sizes: (Custom sizes also available by quotation.)

- 80" x 35" x 6"
- 80" x 36" x 6"
- 80" x 39" x 6"

### Recommended Weight Capacity:

350 lbs (with all foam inserts)

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.

### 100% Customer Satisfaction Guarantee:

The warranty period for this product is 5 years.