

Shape Series Wheelchair Cushions Shape, Visco Shape, and Gel Shape Cushion

Designed for individuals at moderate-to-high risk of skin breakdown. Combines visco and polyurethane foam or gel and polyurethane foam working together, to reduce shear/friction that can potentially lead to pressure ulcers. Lateral supports and an ischial well add comfort and assist in offloading bony prominences.




Includes an inner moisture barrier cover plus a premium quality outer cover.

Key Considerations:

- P** • Best choice for individuals at moderate-to-high risk for pressure sore development
- F** • Recommended for individuals requiring postural support
- T** • Gel option has an added advantage of preventing excessive heat build-up and sweating
- C** • Fully contoured ischial well; leg channeling for thigh positioning
- W** • 250 lbs

ESP Series Wheelchair Cushions ESP or ESP with Viscoelastic Foam

Designed for individuals at low risk of skin breakdown. Proven design using multi-layered latex and polyurethane foams or a combination of viscoelastic and polyurethane foams working together in a low profile cushion.



Lateral supports and an ischial well add comfort and support, ensuring postural stability.

Key Considerations:

- P** • Good choice for individuals at low risk for pressure sore development
- F** • Recommended for individuals requiring modest postural support
- T** • Good choice for prevention of heat build-up and sweating
- C** • Mild leg channeling with waterfall front
- W** • 250 lbs
- O** • Best choice for individuals who foot propel

EMP

Essential Medical Products

www.essentialmedicalproducts.com

B-Gel Wheelchair Cushion Bariatric

Designed for heavier individuals. Combines visco foam and gel working together, reducing shear/friction that can potentially lead to pressure ulcers. Defined ischial well with no leg channeling adds comfort.



Includes an inner moisture-barrier cover plus a premium quality outer cover.

Key Considerations:

- P** • Best choice for individuals at moderate-to-high risk for pressure sore development
- F** • Recommended for individuals requiring postural support
- T** • Gel has an added advantage of preventing excessive heat build-up and sweating
- C** • Fully contoured Ischial well; No leg channeling
- W** • 600 lbs
- O** • Best choice for bariatric applications

Hip Cushion Hip Replacement Recovery

The Arthritis Hip Cushion is a medium density, economical, cushion intended for temporary use after hip replacement surgery.



This light-weight cushion comes complete with a carrying handle to allow the client to carry it with them during their recovery period.

Pin-core Sunmate Wheelchair Cushion

Ideal for Individuals Needing a Softer Cushion

Pillowy soft comfortable layers of latex and viscoelastic foams provide pressure relief.



Ideal for individuals at low-to-moderate risk for pressure sores.


Key Considerations:

- P** • Good choice for individuals at low-to-moderate risk for pressure sore development
- F** • Recommended for individuals requiring a soft cushion
- T** • Good choice for prevention of heat build-up and sweating
- C** • Non-contoured smooth surface
- W** • 250 lbs

Active Series Wheelchair Cushions

Active Firm or Active Plush

An economical choice for individuals requiring a cushion with excellent support qualities, aiding in transfers and repositioning.



Ideal for individuals at low risk for pressure sores.

Key Considerations:

- P** • Good choice for individuals at low risk for pressure sore development
- F** • Available in your choice of plush or firm support
- T** • Good choice for prevention of heat build-up and sweating
- C** • Non-contoured smooth surface
- W** • 250 lbs
- O** • Best choice where follow-up is unlikely. Cushion can be used either side up


P - Pressure, F - Firmness, T - Temperature, C - Contour, W - Weight Capacity, O - Other

604 274 3626

sales@essentialmedicalproducts.com

Anterior Wedge Support

A simple wedge shaped support of firm, closed cell foam material. The wedge can be trimmed to fit over wheelchair seat upholstery up to 20" wide and 20" deep.



Designed to be located under an existing wheelchair cushion. The wedge is 1-1/2" high but can be trimmed to provide less elevation at the front of the seat. Once trimmed the wedge will result in additional anterior wedging of the wheelchair cushion.

Fabricated using a firm, closed cell foam, the finished trimmed wedge may be cleaned in the same way the wheelchair upholstery is cleaned.

Standard Sizes

Depth	Standard Cushions			Extra-Wide Cushions		
	16"	18"	20"	22"	24"	26"
16"	x	x	x			
18"	x	x	x	x	x	x
20"	x	x	x	x	x	x
22"					x	x

Custom Sizes are also available.

Visit www.essentialmedicalproducts.com for further information.