


LTC 3500 Active Care Mattress

Optimal Choice for Ease of Movement

Designed for clients with muscle weakness who find it difficult to move in bed. Combines a top contour cut comfort layer to reduce shear and multi-zoned core enhancing pressure distribution over vulnerable areas.



Key Considerations:

- P** • Functional choice for individuals at risk for pressure sore development
- F** • Recommended for individuals requiring a medium-firm mattress
- H** • Functional choice for prevention of heat build-up and sweating
- T** • Good choice for bedside stability and preventing entrapment
- M** • Best choice for ease of movement on bed surface
- W** • Recommended weight capacity is 350 lbs

LTC 3700 Conformity Mattress

Soft & Cool Memory Foam Mattress

Dual layer mattress designed for individuals looking for envelopment associated with viscoelastic foam while maintaining the ability to move while in bed.



Key Considerations:

- P** • Good choice for individuals at risk for pressure sore development
- F** • Recommended for individuals requiring a medium-soft mattress
- H** • Good choice for prevention of heat build-up and sweating
- T** • Good choice for bedside stability and preventing entrapment
- M** • Viscoelastic top layer may impair bed mobility for some individuals
- W** • Recommended weight capacity is 300 lbs

EMP

Essential Medical Products

www.essentialmedicalproducts.com

LTC 8000 Gel-In Mattress

For Prevention of Heat Build-up and Sweating

Advanced hybrid gel/foam technology with cooling gel infused into the viscoelastic foam top layer which results in a perfect balance of heat regulation and energy absorption. Maximum cooling ability for clients who sweat due to overheating.



Key Considerations:

- P** • Good choice for individuals at risk for pressure sore development
- F** • Recommended for individuals requiring medium firmness
- H** • Best choice for prevention of heat build-up and sweating
- T** • Good choice for bedside stability and preventing entrapment
- M** • Viscoelastic top layer may impair bed mobility for some individuals
- W** • Recommended Weight Capacity is 350 lbs

LTC Bariatric & Ultra Bariatric Mattress

Comfort For the Bariatric Client

- Bariatric mattress: Dual-layer construction designed for heavier clients up to 500 lbs.
- Ultra Bariatric mattress: Tri-layer construction accommodates clients up to 700 lbs.



Key Considerations:


- P** • Good choice for individuals at risk for pressure sore development
- F** • Recommended for individuals requiring a medium-firm mattress
- H** • Functional choice for prevention of heat build-up and sweating
- T** • No transfer borders - reduced stability for unassisted transfers
- M** • Good choice for ease of movement on bed surface
- W** • Best choice for a bariatric mattress

P - Pressure, F - Firmness, H - Heat, T - Transfers, M - Bed Mobility, W - Weight

LTC 4000 Plus Mattress

Pressure Redistribution with Good Mobility

Dual layer mattress designed for individuals looking for a pressure relieving surface that is both supportive and comfortable, combining pressure redistribution with good bed mobility.



Key Considerations:

- P** • Good choice for individuals at risk for pressure sore development
- F** • Recommended for individuals requiring medium firmness
- H** • Good choice for prevention of heat build-up and sweating
- T** • Good choice for bedside stability and preventing entrapment
- M** • Viscoelastic top layer may impair bed mobility for some individuals
- W** • Recommended Weight Capacity is 350 lbs

LTC 4000 Ultra 55 Mattress

Ultra Stable Border to Prevent Entrapment

Triple layer mattress with full height transfer borders ensures stability during transfers and when performing ADL's. Helps prevent entrapment between mattress and bedside rails.



Key Considerations:

- P** • Best choice for individuals at risk for pressure sore development
- F** • Recommended for individuals requiring a firm mattress
- H** • Good choice for prevention of heat build-up and sweating
- T** • Best choice for bedside stability, and for preventing entrapment
- M** • Viscoelastic top layer may impair bed mobility for some individuals
- W** • Recommended Weight Capacity is 400 lbs

P - Pressure, F - Firmness, H - Heat, T - Transfers, M - Bed Mobility, W - Weight

604 274 3626

sales@essentialmedicalproducts.com

Mattress Covers are fluid-proof, non-allergenic, antibacterial, and antimicrobial.

Where **pressure redistribution is a priority**, LTC Emersion (all-way stretch) cover ensures optimal mattress performance.

Need help choosing your mattress?

Visit www.essentialmedicalproducts.com for further information.

Will this mattress be used on a wire spring base platform?

You need a Mattress Protector (MBP 500)

Is patient in-bed safety a concern?

Look for LTC Side Rail Pads, Headboard and Footboard Pads

Are patient falls a concern?

Look for LTC Bedside Security Pads

