

Gel-Shape Wheelchair Cushion

The **Gel-Shape Wheelchair Cushion** incorporates specialized Gel technology which allows this cushion to conform easily to the body for freedom of movement that is far superior to common padding products. Its unique design prevents friction between the padding material and human skin to help prevent damage to tissue and nerves that causes skin breakdown. Lateral supports and a central pommel help to evenly distribute weight for increased comfort and stability.



- Incorporates specialized gel technology
- Prevents damage to tissue and nerves
- Minimizes friction
- Evenly distributes weight

The incontinent mesh-stretch cover is non-allergenic, bacteria resistant and wipes clean to help prevent cross infections. A non-skid bottom prevents cushion travel.

A specially priced, additional cover is available at time of purchase.

Standard Sizes: (Custom sizes also available by quotation.)

- 16" x 16"
- 16" x 18"
- 16" x 20"
- 18" x 16"
- 18" x 18"
- 18" x 20"

Recommended Weight Capacity:

300 lbs

Height:

3 3/4" contoured surface

100% Customer Satisfaction Guarantee:

For over 35 years, Essential Medical Products has been committed to using the highest quality, most environmentally friendly materials available, in all of its products. Our goal is to enhance the quality of care and guarantee each customer's satisfaction.

Adding Comfort to Your Care