

LTC 3500 Active Care Mattress

Designed for clients with muscle weakness who find it difficult to move in bed. Combines a top contour cut comfort layer to reduce shear and multi-zoned core enhancing pressure distribution over vulnerable areas.




Patient Recommendations:

- Functional choice for individuals at risk for pressure sore development
- Recommended for individuals requiring a medium-firm mattress
- Functional choice for prevention of heat build-up and sweating
- Good choice for bedside stability, preventing entrapment
- Best choice for ease of movement on bed surface
- Recommended weight capacity is 350 lbs

LTC 3700 Conformity Mattress Soft & Cool Memory Foam Mattress

Dual layer mattress designed for individuals looking for envelopment associated with viscoelastic foam while maintaining the ability to move while in bed.



Patient Recommendations:

- Good choice for individuals at risk for pressure sore development
- Recommended for individuals requiring a medium-soft mattress
- Good choice for prevention of heat build-up and sweating
- Good choice for bedside stability, preventing entrapment
- Viscoelastic top layer may impair bed mobility for some individuals
- Recommended weight capacity is 300 lbs

EMP


Essential Medical Products

www.essentialmedicalproducts.com

LTC 1500 Mattress

Basic Care Mattress

Designed for Long Term Care general purpose applications. Most economical LTC mattress available.




Patient Recommendations:

- Not recommended for individuals at risk for pressure sore development
- Recommended for individuals requiring a medium-soft mattress
- Functional choice for prevention of heat build-up and sweating
- No transfer borders - reduced stability for unassisted transfers
- Good choice for ease of movement on bed surface
- Most economical choice for general applications
- Recommended Weight Capacity is 250 lbs

LTC Bariatric & Ultra Bariatric Mattress Comfort For the Bariatric Client

- Bariatric mattress: Dual-layer construction designed for heavier clients up to 500 lbs.
- Ultra Bariatric mattress: Tri-layer construction accommodates clients up to 700 lbs.



Patient Recommendations:

- Good choice for individuals at risk for pressure sore development
- Recommended for individuals requiring a medium-firm mattress
- Functional choice for prevention of heat build-up and sweating
- No transfer borders - reduced stability for unassisted transfers
- Good choice for ease of movement on bed surface
- Best choice for bariatric individuals
- Recommended Weight Capacity is:
Bariatric Mattress - 500 lb / Ultra Bariatric Mattress - 700

LTC 4000 Ultra 55 Mattress

Ultra Stable Firm Memory Foam

Firm and deep viscoelastic foam with full height transfer borders ensure stability during transfers and when performing ADL's. Helps prevent entrapment between mattress and bedside rails.



Patient Recommendations:

- Best choice for individuals at risk for pressure sore development
- Recommended for individuals requiring a firm mattress
- Good choice for prevention of heat build-up and sweating
- Best choice for bedside stability, preventing entrapment
- Thick viscoelastic top layer may impair bed mobility for some individuals
- Good choice for heavier individuals or individuals preferring a firm mattress
- Recommended Weight Capacity is 450 lbs

LTC 8000 Gel-In Mattress

For Heat Regulation and Energy Absorption

Advanced hybrid gel/foam technology with cooling gel infused into the viscoelastic foam top layer which results in a perfect balance of heat regulation and energy absorption. Maximum cooling ability for clients who sweat due to overheating.



Patient Recommendations:

- Best choice for individuals at risk for pressure sore development
- Recommended for individuals requiring medium firmness
- Best choice for prevention of heat build-up and sweating
- Good choice for bedside stability and preventing entrapment
- Viscoelastic top layer may impair bed mobility for some individuals
- Recommended Weight Capacity is 350 lbs

604 274 3626

sales@essentialmedicalproducts.com

Need help choosing your mattress?

Visit www.essentialmedicalproducts.com for further information.

Will this mattress be used on a wire spring base platform?

You need a Mattress Protector (MBP 500)

Is patient in-bed safety a concern?

Look for LTC Side Rail Pads, Headboard and Footboard Pads

Are patient falls a concern?

Look for LTC Bedside Security Pads

Patient positioning in bed?

Look for LTC Body Positioning Wedges and Heel Lift Wedges.

Need other bed safety products?

Look for LTC Care Aid Bed Positioner.

