



Mattress Covers – Choose which cover material best meets your patients needs

In many cases you may now choose the type of cover material for your mattress. All EMP mattress covers are fluid-proof, non-allergenic, antibacterial, antimicrobial and tested to ensure long-life in the hospital and home care environments.

- **Lux NS** is the cover material supplied with all EMP mattresses for the past 15+ years. This material enjoys an excellent reputation as the most durable cover available for continuing care applications
- **Emersion** is a new stretch material developed at EMP to provide improved pressure redistribution characteristics, allowing both urethane and viscoelastic foams to perform optimally. This material is a true 4-way stretch fabric that is PVC free, fire rated, and incorporates ultrasonic welded seam technology, which is better from both infection and pest control standpoints.

The table below identifies which mattresses can be ordered with **Lux NS** (non-stretch) and/or **Emersion** (4-way stretch) cover material.

Description	LTC 1500	LTC 2500	LTC 3500	LTC 3700	LTC 4000S	LTC 4000U55
Mattress Name	Basic Care	Basic Comfort	Active Care	Conformity	Plus	Ultra 55
Standard Cover	Lux NS	Lux NS	Emersion	Emersion	Emersion	Emersion
	Non-Stretch	Non-Stretch	Stretch	Stretch	Stretch	Stretch
Standard Cover Product ID	Std Code	Std Code	Std Code	Std Code	Std Code	Std Code
Alternate Cover	Not Available	Not Available	Lux NS	Lux NS	Lux NS	Lux NS
			Non-Stretch	Non-Stretch	Non-Stretch	Non-Stretch
Alternate Cover Product ID (add when ordering)	N/A	N/A	Suffix "NS"	Suffix "NS"	Suffix "NS"	Suffix "NS"

Description	LTC 8000	LTC 9000	LTC 9100	9100 Inserts	LTC Bariatric	LTC Ultra Bariatric
Mattress Name	Gel-In	Progressive Care	Intensive Care		Bariatric	Ultra Bariatric
Standard Cover	Emersion	Emersion	Emersion	Emersion	Lux NS	Lux NS
	Stretch	Stretch	Stretch	Stretch	Non-Stretch	Non-Stretch
Standard Cover Product ID	Std Code	Std Code	Std Code	Std Code	Std Code	Std Code
Alternate Cover	Lux NS	Lux NS	Lux NS	Lux NS	Emersion	Emersion
	Non-Stretch	Non-Stretch	Non-Stretch	Non-Stretch	Stretch	Stretch
Alternate Cover Product ID (add when ordering)	Suffix "NS"	Suffix "NS"	Suffix "NS"	Suffix "NS"	Suffix "EM"	Suffix "EM"

Please contact us for help determining the best covers for your facility:

Phone: 604 274 3626 / 1 800 344 4437

sales@essentialmedicalproducts.com

www.essentialmedicalproducts.com