



EMP

Essential Medical Products

Adding Comfort to Your Care

Cushions

Clinician's choice when quality and performance matter most!




Shape Series Wheelchair Cushions Shape, Visco Shape, and Gel Shape Cushion

Designed for individuals at moderate-to-high risk of skin breakdown. Combines polyurethane foam, viscofoam or gel working together, to reduce shear/friction that can potentially lead to pressure ulcers. Lateral supports and an ischial well add comfort and assist in offloading bony prominences.

Includes an inner moisture barrier cover plus a premium quality incontinent outer cover.

Key Considerations:

- P** • Best choice for individuals at moderate risk for pressure sore development
- F** • Recommended for individuals requiring postural support
- T** • Gel option has an added advantage of preventing excessive heat build-up and sweating
- C** • Fully contoured ischial well; leg channeling for thigh positioning
- W** • 250 lbs
- O** • Best choice for maximum support



P - Pressure, F - Firmness, T - Temperature, C - Contour, W - Weight Capacity, O - Other




ESP Series Wheelchair Cushions ESP/ESP Deluxe with Viscoelastic Foam

Designed for individuals at low risk of skin breakdown. Proven design using multi-layered latex and polyurethane foams or a combination of viscoelastic and polyurethane foams working together in a low profile cushion.

Lateral supports and an ischial well add comfort and support, ensuring postural stability. Deluxe version includes inner moisture barrier. Either version includes incontinent outer cover.

Key Considerations:

- P** • Good choice for individuals at low risk for pressure sore development
- F** • Recommended for individuals requiring modest postural support
- T** • Good choice for prevention of heat build-up and sweating
- C** • Mild leg channeling with waterfall front
- W** • 250 lbs
- O** • Best choice for individuals needing a low seat height and/or who foot propel



P - Pressure, F - Firmness, T - Temperature, C - Contour, W - Weight Capacity, O - Other

www.essentialmedicalproducts.com



604 274 3626

sales@essentialmedicalproducts.com




Pincore Viscoelastic Wheelchair Cushions Ideal for Individuals Needing a Softer Cushion

Pillowy soft comfortable layers of latex and viscoelastic foams provide pressure redistribution.

Ideal for individuals at low-to-moderate risk for pressure sores. Incontinent cover.

Key Considerations:

- P** • Good choice for individuals at low-to-moderate risk for pressure sore development
- F** • Recommended for individuals requiring a soft cushion
- T** • Good choice for prevention of heat build-up and sweating
- C** • Non-contoured smooth surface
- W** • 250 lbs
- O** • Extremely comfortable all-purpose cushion



Active Series Wheelchair Cushions Active Firm or Active Plush

An economical choice for individuals requiring a cushion with excellent support qualities, aiding in transfers and repositioning.

Ideal for individuals at low risk for pressure sores. Includes incontinent cover.

Key Considerations:

- P** • Good choice for individuals at low risk for pressure sore development
- F** • Available in your choice of plush or firm support
- T** • Good choice for prevention of heat build-up and sweating
- C** • Non-contoured smooth surface
- W** • 250 lbs
- O** • Best choice where clinical follow-up is less frequent
Cushion can be used either side up



Standard Sizes:

Depth	Standard Cushions			Shape XL Cushions		
	16"	18"	20"	22"	24"	26"
16"	x	x				
18"	x	x	x	x	x	x
20"	x	x	x	x	x	x
22"					x	x

Additional Support - Anterior Wedge Support

A simple wedge shaped support of firm, closed cell foam material. The wedge can be trimmed to fit over wheelchair seat upholstery up to 20" wide and 20" deep. Designed to be located under an existing wheelchair cushion. The wedge is 1-1/2" high but can be trimmed to provide less elevation at the front of the seat. Once trimmed the wedge will result in additional anterior wedging of the wheelchair cushion.

Fabricated using a firm, closed cell foam, the finished trimmed wedge may be cleaned in the same way the wheelchair upholstery is cleaned.



EMP
Essential Medical Products



Enhanced
Comfort
and Support



Longlasting
Quality
Foam Cushions

Use this guide to help choose the best option for individual patients.

Essential Medical Products offers a comprehensive selection of standard and custom medical products for bedding, rehabilitation therapy and wheelchair applications.

Contact us for superior medical cushioning products:

Phone: 604 274 3626 / 1 800 344 4437 Fax: 604 274 3627 / 1 800 438 3381
sales@essentialmedicalproducts.com www.essentialmedicalproducts.com

