


## Conform<sup>®</sup> Air Insert

The ideal support surface environment when pressure and/or moisture is of concern.

Designed to redistribute applied pressure evenly across a surface area, Conform<sup>®</sup> Inserts ensure normally high-risk pressure points are minimized or eliminated. Each air insert includes a grid pattern of individual bulb-shaped supports with space between each, ensuring air can flow, allowing moisture to move away from the skin surface.




### Patient Considerations for this Insert:

Excellent pressure relief - encouraging capillary blood flow. Can be used from the stage of prophylaxis up to pressure ulcer (Grade IV – Highest Risk). Waterlow score - these products can be used at a score of 20 and above.

### CONFORM<sup>®</sup> Med. Air Inserts are made of allergen-free natural rubber.

Available in two sizes, ensuring a perfect fit in a variety of mattress applications. EMP recommends inserts are best used with EMP Progressive Care and Intensive Care mattresses. Visit [www.essentialmedicalproducts.com](http://www.essentialmedicalproducts.com) for more information.



### Recommended Weight Capacity:

350 lbs  
Bariatric applications are permitted with a shorter warranty period.

### 100% Customer Satisfaction Guarantee:

The warranty period is 30 months.  
Bariatric applications warranty period will be 12 months.

Product Code	Size
CAS-30W	30" x 19" x 2-3/4"
CAS-34W	34" x 19" x 2-3/4"



Connector Strip used to join two or more inserts

This information provides general guidance. EMP recommends specialized sleep surfaces be selected with the help of a qualified medical professional.