


LTC 8000 Gel-In Mattress

For Heat Regulation and Energy Absorption

Advanced hybrid gel/foam technology with cooling gel infused into the viscoelastic foam top layer which results in a perfect balance of heat regulation and energy absorption. Maximum cooling ability for clients who sweat due to overheating.



Patient Considerations for this Mattress:

- | | |
|------------------------------|---|
| Pressure Distribution | › Best choice for individuals at risk for pressure sore development |
| Firmness | › Recommended for individuals requiring medium firmness |
| Temperature | › Best choice for prevention of heat build-up and sweating |
| Transfers | › Best choice for bedside stability and preventing entrapment |
| Bed Mobility | › Viscoelastic top layer may impair bed mobility for some individuals |

Specifications:

- | | |
|-------------------------|---|
| Top Foam Layer | › 2" gel infused viscoelastic foam comfort layer |
| Core Foam | › Multi-zone medium-firm core |
| Transfer Borders | › 6" high exterior foam adds firmness along edges |
| Cover | › Fluid-proof, non-allergenic, antibacterial, antimicrobial cover
Also available with organic cotton cover (not fluid-proof) |
| Misc | › Designed for use on all hospital beds |

Standard Sizes: (Custom sizes also available by quotation.)

- | | |
|------------------|------------------|
| • 80" x 35" x 6" | • 80" x 48" x 6" |
| • 80" x 36" x 6" | • 80" x 54" x 6" |
| • 80" x 39" x 6" | • 80" x 59" x 6" |

Recommended Weight Capacity:

350 lbs

100% Customer Satisfaction Guarantee:

The warranty period is five years; two years full plus three years pro-rated.

This product is intended for use by persons with a permanent disability. EMP recommends all EMP products be selected with the help of a qualified medical professional.