


# EMP

Essential Medical Products



## Conform<sup>©</sup> Air Segment

The ideal support surface environment when pressure injury and/or moisture is of concern.

Designed to redistribute applied pressure evenly across a surface area, Conform<sup>©</sup> Segments ensure normally high-risk pressure points are minimized or eliminated. Each segment includes a grid pattern of individual bulb-shaped supports with space between each, ensuring air can flow, allowing moisture to move away from the skin surface.

### Patient Considerations for this Insert:

Excellent pressure relief - encouraging capillary blood flow. Can be used from the stage of prophylaxis up to pressure ulcer (Grade IV – Highest Risk). Waterlow score - these products can be used at a score of 20 and above.

**CONFORM<sup>®</sup> Med. Air Segments are made of allergen-free natural rubber.**

### Recommended Weight Capacity:

350 lbs

Bariatric applications are permitted with a shorter warranty period.

### 100% Customer Satisfaction Guarantee:

The warranty period is 30 months.

Bariatric applications warranty period will be 12 months.



Connector Strip used to join two or more segments

Product Code	Size
CAS-30	30" x 19" x 2-3/4"
CAS-34	34" x 19" x 2-3/4"

This product is intended for use by persons with a permanent disability. EMP recommends all EMP products be selected with the help of a qualified medical professional.