

Premium Care Mattresses

Where Do I Use an EMP Premium Care Mattress?

- Premium Care mattresses are 'passive support surfaces', ideal for use at home, in care, or in hospital
- Designed for individuals at high risk of developing a pressure injury, or with a Stage I injury

What options can I select?

- Transfer borders (when included) provide firm support when getting in or out of bed
- Choose from a variety of high-quality inserts that provide pressure relief, improved airflow and patient support

Why Prescribe an EMP mattress?

- Higher density foam results in durability of feel each and every night
- CertiPUR certification ensures the mattress is free of toxic materials. Mattress contains Natural Oil Polyols.
- Proven reputation for quality
- Superior mattress cover outlasts the competition
- Cover easy to clean using sanitary wipes
- 3" topper thickness for optimal immersion in areas not filled using mattress inserts
- 3" insert thickness for optimal immersion
- Available in 35", 36", 39", 42" 48" and 60" wide. 36" wide always in-stock for fast delivery

EMP

Essential Medical Products


Adding Comfort to Your Care

LTC 9000 Progressive Care Mattress
Adjustable Pressure Redistribution

Inserts are UNCOVERED. The entire mattress is enveloped in 4 way stretch (Emersion) cover.



Available with 3" or 6" tall transfer borders



Insert Choices:


- Uncovered conform inserts
- Uncovered foam inserts

LTC 9100 Intensive Care Mattress
Pressure Redistribution Focused on Specific High Risk Areas

Foam inserts are COVERED in 4 way stretch (Emersion) material. There are options for friction fit pads; snaps and flaps, or Velcro flaps.



Available with 3" or 6" tall transfer borders



Insert Choices:


- Uncovered or covered conform inserts
- Covered foam inserts

LTC 9200 T-Style Configurable Mattress
Pressure Re-distribution Focused on Coccyx, Sacrum and Heels

Foam inserts are COVERED in 4 way stretch (Emersion) material. Velcro flaps secure foam inserts to the mattress base.



Transfer Rails located under Viscoelastic Top Layer



Bird's-eye View


Insert Choices:

- Uncovered or covered conform inserts
- Covered foam inserts


Choose from the following high-quality Inserts:

- Conform Air Inserts - provide the best passive pressure redistribution as well as air flow at the skin surface
- Latex foam Inserts - naturally conforms to the body's contours for optimal pressure redistribution (as in the LTC 7000 mattress)
- Gel-in memory foam Inserts - ideal for individuals needing a surface having pressure relieving and cooling properties (as in the LTC 4000 Plus mattress)
- Polyurethane foam Inserts - used in situations where you want insert 'slots' to be available for future use when needed (as in the LTC 3500 mattress)


Conform Air Inserts (covered or uncovered)




















Covered Foam Inserts



Uncovered Foam Inserts



Premium Care Inserts At-A-Glance

|  EMP Essential Medical Products | Polyurethane Foam Inserts | Gel-In Memory Foam Inserts | Latex Foam Inserts | Conform Air Inserts |
|--|---|--|---|---|
| Pressure Redistribution |  |  |  |  |
| Firmness |  |  |  |  |
| Temperature |  |  |  |  |
| Mobility |  |  |  |  |
| Weight Capacity | 350 lbs | 350 lbs | 350 lbs | 350 lbs |

TEMPERATURE (Gauge 1=warm and 10=cool) Illustrates the foam's ability to draw heat away from the body, or to hold the body's heat close to its source. The higher the gauge reads, the cooler your body's environment will stay throughout the night.

FIRMNESS (Gauge 1=soft and 10=firm) This gauge is a direct reading of how firm, or soft, the sleep surface will be. The higher the reading on the gauge, the firmer the sleep surface.

This product is intended for use by persons with a permanent disability. EMP recommends all EMP products be selected with the help of a qualified medical professional.

PRESSURE REDISTRIBUTION (Gauge 1=less and 10=more) This gauge illustrates how much of the patient's total body weight is redistributed through the foam, thus alleviating pressure points. The higher the reading on this gauge indicates that the weight placed on this foam is spread out over a larger area, thus decreasing the overall pressure felt by the body.

MOBILITY (Gauge 1=less and 10=more) The mobility rating represents the foam's ability to "push-back" against the weight load that is on its surface. The higher the rating on this gauge, the easier it is for a person to move on this foam type.

Contact us for superior medical cushioning products:

Phone: 604 952 7005 / 1 800 344 4437 Fax: 604 952 7035 / 1800 438 3381
cs@essentialmedicalproducts.com www.essentialmedicalproducts.com

Essential Medical Products Ltd. is a subsidiary of BFF Foam and is proudly 100% Canadian.

REV 06/2022

